

Name:

Date:

Class:

### Field of View Student Activity Sheet

A. Begin using low power. Put a clear plastic ruler across the hole in the stage. Look through the microscope and find the ruler's millimeter (mm) marks. Count the number of mm that fit *across* the field of view. (Note: *Field of view* is a term that refers to "everything you see" when looking through the eye piece of a microscope.)

- Measure: \_\_\_\_\_ mm (all measurements should go to three digits)
- Now multiply by 1000 to convert mm to micrometers ( $\mu\text{m}$ ). \_\_\_\_\_  $\mu\text{m}$ .

This is the size of your field of view.

B. Move to medium power.

- Count the number of mm that fit across the field of view. Be precise!  
Measure \_\_\_\_\_ mm
- Now multiply by 1000 to convert mm to micrometers ( $\mu\text{m}$ ): \_\_\_\_\_  $\mu\text{m}$ .

What can you say about the relationship between the *degree* of magnification and the *size* of your field of view?

C. Move to high power.

- Notice that field of view is now *less* than a millimeter across. In order to calculate the size of your field of view, you need to use this formula:

$$\frac{\text{size of the field on low power, in } \mu\text{m}}{\text{Total magnification on low power}} \times \text{Total magnification on high power} = \text{size of the field of view on high power } \mu\text{m}.$$

Total magnification on high power

size of the field of view on high power

Questions:

1. Using whatever magnification power is most convenient, estimate the width of a newspaper letter "e", in  $\mu\text{m}$ , by comparing it to the overall size of the field of view.

Answer: \_\_\_\_\_  $\mu\text{m}$

2. Why do you think it is important to know the size of the field of view?